## EDTA

EDTA | EDTA  
Ref. 35.001.00

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### ENGLISH

#### INTENDED USE

Anticoagulant solution for hematology tests and plasma measurements.

#### PRINCIPLE OF THE METHOD

The EDTA (ethylenediamine tetracetic acid) is used to obtain total blood and plasma. The EDTA acts chelating the ionic calcium leading to formation of EDTA-Calcium complex. Thus cascade effect of coagulation (calcium dependent) stops, blocking the conversion of Prothrombin into Thrombin and the interaction of the Thrombin with Fibrinogen in the Fibrin (clot).

#### REAGENT COMPOSITION

<b>R 1</b>	EDTA sodium 300 mmol/L, met hyl orange 22 mg/L
------------	--

#### STORAGE AND STABILITY

- The reagents are stable up to the date stated on the label.
- Do not use reagents over the expiration date.
- Store at 15 to 30 °C.

#### ADDITIONAL EQUIPMENT

- Assay tubes and the usual material for blood collection.

### WARNINGS AND PRECAUTIONS

- The kit is intended only for in vitro diagnostic use.
- The EDTA Anticoagulant must not be used in ions or prothrombin time assays.
- The safety procedures for handling biological and immunochemical materials must be observed.
- Wear the individual protection equipment according to Good Laboratory Practice.
- Disposal of all waste material should be in accordance with local guidelines.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling for 55-35-3214-4646.
- In case of ingestion seek for medical advice.
- The reagents must be kept out of the specified temperatures only the time necessary for the realization of the tests.
- Do not combine reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign of contamination as turbidity, flocculation and/or color change.

### SAMPLE – PREPARATION AND STABILITY

Blood

#### TECHNICAL PROCEDURE

##### A) REAGENT PREPARATION

Ready to use.

##### B) PROCEDURE

- Use 50µL (one drop) of anticoagulant for each 5 mL of blood.
- Homogenize gently by inversion.

Note: The anticoagulant/blood ratio must be observed in order to avoid hemolysis.

#### TECHNICAL LIMITATIONS

Hemolysis may occur in cases of anticoagulant excess, inappropriate homogenization of the sample and cleaning of lab material and/or when difficulties in the blood collection occur.

#### QUALITY CONTROL

Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Asq, it must be kept a defined system for verification of the analytical variability that occurs in any measuring system. The lab must participate of external quality control programs.

#### OBSERVATIONS

The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.

2. The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this

product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotecnica.com.br](mailto:sac@biotecnica.com.br)

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