## FLUORETO

Fuori de | Fluoruro  
Ref. 35.003.00

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### ENGLISH

#### INTENDED USE

Anticoagulant solution for glycolysis inhibition.

#### PRINCIPLE OF THE METHOD

The EDTA in the anticoagulant acts blocking the ionic calcium leading for forming an EDTA-calcium complex inhibiting the coagulation and the fluoride inhibits the glycolysis.

#### REAGENT COMPOSITION

<b>R 1</b>	Potassium fluoride 0,25 mmol/L; EDTA 180 mmol/L and light green 50 ng/L
------------	---

#### STORAGE AND STABILITY

- The reagents are stable up to the date stated on the label.
- Do not use reagents over the expiration date.
- Store at 15 to 30 °C and protect from light.

#### ADDITIONAL EQUIPMENT

- Assay tubes and the usual material for blood collection.

#### WARNINGS AND PRECAUTIONS

- The kit is intended only for in vitro diagnostic use.

- The safety procedures for handling biological and immunochemical materials must be observed.
- Wear the individual protection equipment according to Good Laboratory Practice.
- Disposal of all waste material should be in accordance with local guidelines.
- In case of ingestion seek for medical advice.
- The reagents must be kept out of the specified temperatures only when necessary for the realization of the tests.
- Do not combine reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign of contamination as turbidity, flocculation and/or color change.

#### SAMPLE – PREPARATION AND STABILITY

Blood.

#### TECHNICAL PROCEDURE

##### A) REAGENT PREPARATION

Ready to use.

##### B) PROCEDURE

- Use 50 µL (one drop) of anticoagulant for each 3 mL of blood.
- Homogenize gently by inversion.
- Centrifuge at 3.000 rpm for 5 minutes.
- Transfer the supernatant to another tube and keep it refrigerated up to use.

Note: The anticoagulant/blood ratio must be observed to avoid hemolysis.

#### TECHNICAL UTILIZATIONS

Hemolysis may occur in cases of anticoagulant excess, inappropriate homogenization of the sample and cleaning of lab material and/or when difficulties in the blood collection occur.

#### QUALITY CONTROL

Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Asq it must be kept a defined system for verification of the analytical variability that occurs in any measuring system.

The lab must participate of external quality control programs.

#### OBSERVATIONS

The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.

The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this

product or this procedure, contact us calling +55 35 3214 4646, your local distributor or send an e-mail: [sac@biotecnica.com.br](mailto:sac@biotecnica.com.br)

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